



### Cleanroom Gloves

#### TRIonic® Z-Grip Gloves

Ideal for use in mixing, handling, and transporting chemicals. They feature a "Z" pattern to enhance grip when handling wet items. Made of natural rubber, nitrile, and neoprene. Bag of 12 pairs. Individually packaged per pair.

**Approval:** CE

**Compliance:** Class 10 compatible; exceeds ASTM 3577-91, 1.5% AQL; 21CFR 170-199; finished, inspected, and packaged in certified M 3.5 environment.



Catalog number	Size	Length	Thickness	Price per bg
<a href="#">A-40214-00</a>	6-6½	14"L	20 mil	
<a href="#">A-40214-01</a>	7-7½			
<a href="#">A-40214-02</a>	8-8½			
<a href="#">A-40214-03</a>	9-9½	14"L	20 mil	
<a href="#">A-40214-04</a>	10-10½			
<a href="#">A-40214-05</a>	11			

#### TRIonic® Diamond Grip Gloves

Gloves provide superior performance when using aggressive chemicals. The raised diamond grip lets you safely handle wet items. Made of natural rubber, nitrile, and neoprene. Bag of 12 pairs, individually packaged per pair.

**Compliance:** Class 10 compatible; meets or exceeds ASTM 3577-91, 1.5% AQL; 21 CFR 170-199; finished, inspected, and packaged in certified M 3.5 environment



Catalog number	Size	Length	Thickness	Price/bg
<a href="#">A-40214-30</a>	7	14"L	20 mil	
<a href="#">A-40214-31</a>	8			
<a href="#">A-40214-32</a>	9			
<a href="#">A-40214-33</a>	10			
<a href="#">A-40214-34</a>	11			



### StanSolv® Gloves

Nitrile formulation gloves have excellent cut, puncture, and abrasion resistance. Use for handling organic solvents and chemicals and chemical mixing. The "Z" pattern enhances grip when handling wet items. Bag of 12 pairs, individually packaged per pair.

**Approval:** CE

**Compliance:** Class 10 compatible; ISO 9002; meets or exceeds ASTM 3577-91, 1.5% AQL; 21 CFR 170-199; finished, inspected, and packaged in certified M 3.5 environment



Catalog number	Size	Length	Thickness	Price/bg
<a href="#">A-40214-20</a>	6-6½	13"L	11 mil	
<a href="#">A-40214-21</a>	7-7½			
<a href="#">A-40214-22</a>	8-8½			
<a href="#">A-40214-23</a>	9-9½	13"L	11 mil	
<a href="#">A-40214-24</a>	10-10½			
<a href="#">A-40214-25</a>	11			

### Nitrilite® Gloves

Lightweight yet rugged gloves give you increased comfort and dexterity without losing durability. They have a puncture resistance three times greater than vinyl gloves. Nonslip finish; textured tips. Bag of 50 pairs.

**Approval:** CE **Compliance:** Class 100

Catalog number	Size	Length	Thickness	Price/bg
<a href="#">A-33671-00</a>	Small	12"L	3 mil	
<a href="#">A-33671-10</a>	Medium			
<a href="#">A-33671-20</a>	Large			
<a href="#">A-33671-30</a>	X-large			



### Nytek® Gloves

Gloves are made of solvent-resistant nitrile rubber—permeation tested to ASTM standards. Individual pairs are heat-sealed in a polybag to protect gloves during transit and storage. Extended length and roll-up antidrip cuffs protect forearms. Nonslip finish. Pack of 12 pairs.

**Approval:** CE **Compliance:** Class 100

Catalog number	Size	Length	Thickness	Price/pk
<a href="#">A-09717-74</a>	6-6½	13¾"L	15 mil	
<a href="#">A-09717-76</a>	7-7½			
<a href="#">A-09717-78</a>	8-8½			
<a href="#">A-09717-80</a>	9-9½	13¾"L	15 mil	
<a href="#">A-09717-82</a>	10			
<a href="#">A-09717-84</a>	11			



### 100% Nylon Knit Glove Liners

Experience the softness and comfort of seamless, fine-gauge nylon knit. Liners prevent irritation during extended wear of protective gloves. Dust-free nylon is also ideal for cleanrooms. Order full-fingered liners for total protection or half-fingered liners for greater tactile sensitivity. Wash and sterilize for repeated use. Pack of 12 pair.



Catalog number	Size	Price per pack	
		1 pk	10 pk
<b>Full-finger glove liners</b>			
<a href="#">BH-86382-00</a>	Small Medium Large X-large		
<a href="#">BH-86382-01</a>			
<a href="#">BH-86382-02</a>			
<a href="#">BH-86382-03</a>			
<b>Half-finger glove liners</b>			
<a href="#">BH-86382-10</a>	Small Medium Large X-large		
<a href="#">BH-86382-11</a>			
<a href="#">BH-86382-12</a>			
<a href="#">BH-86382-13</a>			