

**buy more,
save more**

6-29 packs **save 10%** 30-119 packs **save 15%** 120+ packs **save 20%**

Mix and match packs of any items from the "blue" tables to save!



Cutter/dispenser
06720-34



Parafilm® Wrap and Dispensers

Colorless, moistureproof plastic film molds to most surfaces. Seal beakers and flasks; wider size covers trays and basins. Film 06720-25 comes in a mailing tube; films 06720-40 and -50 come in a paper dispenser box.

| Catalog number | Width | Length | Price/rl |
|----------------|-------|-----------|----------|
| K-06720-25 | 20" | 50 ft/rl | |
| K-06720-40 | 4" | 125 ft/rl | |
| K-06720-50 | 2" | 250 ft/rl | |

K-06720-34 Parafilm dispenser clear acrylic



06726-00

Plastic Film, PVC

Clear, self-sealing wrapping film. Packaged in a dispenser box with serrated cutting strip.

| Catalog number | Dimensions | Qty/pk | Price/1 pk | Save 10% 6-29 pks | Save 15% 30-119 pks | Save 20% 120+ pks |
|----------------|---------------|--------|------------|----------------------|------------------------|----------------------|
| K-06276-00 | 18" x 2000 ft | 1 | | | | |

GO to page(s) 161 to 168

For a complete line of our beakers, see pages 161 to 168.



06275-00

Aluminum Foil

All-purpose foil used for freezer storage, wrapping, sealing, autoclaving, lining, and weighing or transferring samples. Packaged in a dispenser box with serrated cutting strip.

| Catalog number | Thickness | Dimensions | Qty/pk | Price/1 pk | Save 10% 6-29 pks | Save 15% 30-119 pks | Save 20% 120+ pks |
|----------------|------------|--------------|--------|------------|----------------------|------------------------|----------------------|
| K-06275-00 | Standard | 18" x 500 ft | 1 | | | | |
| K-06275-50 | Heavy-duty | 18" x 500 ft | 1 | | | | |



New

06275-02

Polyfoil™ Square Wraps

Polyester-laminated foil squares provide a chemically inert barrier to cover vessels for autoclaving, incubating, and storage. Ideal for use with organic solvents and caustic solutions. Easy to apply and remove, and fold easily without kinking. Resistant to punctures, tears, and corrosion. **Autoclavable.**

| Catalog number | Material | Dimensions | Qty/pk | Price/1 pk | Save 10% 6-29 pks | Save 15% 30-119 pks | Save 20% 120+ pks |
|----------------|----------------|-----------------|--------|------------|----------------------|------------------------|----------------------|
| K-06275-02 | Foil/polyester | 4" x 4" | 100 | | | | |
| K-06275-04 | | 5 1/8" x 5 1/8" | | | | | |
| K-06275-06 | | 8" x 8" | | | | | |



34502-90



34503-90

PYREX® Glassware Starter Sets

These glassware starter sets contain one each of the most popular sizes. The 1000 Griffin beaker assortment includes 50-mL, 100-mL, 250-mL, 600-mL, and 1-L beakers. The 4980 Erlenmeyer flask assortment includes 50-mL, 125-mL, 250-mL, 500-mL, and 1-L flasks.

| Catalog number | Description | Price |
|----------------|-----------------------------------|-------|
| EG-34502-90 | 1000 Griffin beaker starter set | |
| EG-34503-90 | 4980 Erlenmeyer flask starter set | |