

## A+ Aloe™ Nitrile Gloves

### Formulated with skin-protecting aloe

Powder-free gloves feature a unique aloe coating on the interior glove surface for a comfortable, soothing feel against bare hands. Aloe prevents the drying and chapping that plagues frequent glove wearers.

Gloves are latex-free and contain no natural rubber proteins. Synthetic formulation offers good chemical and puncture resistance. Beaded cuff ensures easy donning. Ambidextrous.



86480-10

### Specifications & Ordering Information

| Catalog number              | Size    | Length | Thickness | Qty/bx | Price/box |
|-----------------------------|---------|--------|-----------|--------|-----------|
| <a href="#">KH-86480-10</a> | Small   | 9½"    | 5 mil     | 100    |           |
| <a href="#">KH-86480-12</a> | Medium  |        |           |        |           |
| <a href="#">KH-86480-14</a> | Large   |        |           |        |           |
| <a href="#">KH-86480-16</a> | X-large |        |           |        |           |

[KH-86402-00](#) Wire glove dispenser, single

**GO to page(s) 1479-1484**

We offer a complete supply of lab coats, coveralls, aprons, sleeves, and boot covers. See our "Safety-Apparel" section.

## FreeForm™ SE Nitrile Gloves

### Patented "super soft" formulation feels like latex

Powder-free, 100% nitrile gloves offer superior comfort and elasticity in a synthetic glove. Unique formulation conforms to the contours of your hands as you move—minimizing unnecessary hand stress and enabling a full range-of-motion.

Lightly textured fingers ensure a secure grip in wet and dry conditions. Interior polymer coating assists in quick and easy donning.

Compliance: ASTM D 6319.



86445-01

### Specifications & Ordering Information

| Catalog number              | Size    | Length | Thickness | Qty/bx | Price/box |
|-----------------------------|---------|--------|-----------|--------|-----------|
| <a href="#">KH-86445-00</a> | Small   | 9½"    | 6 mil     | 100    |           |
| <a href="#">KH-86445-01</a> | Medium  |        |           |        |           |
| <a href="#">KH-86445-02</a> | Large   |        |           |        |           |
| <a href="#">KH-86445-03</a> | X-large |        |           |        |           |

## Touch N Tuff™ Nitrile Gloves

### The standard in comfort, durability, and cleanliness

Premium quality gloves are thin enough to allow excellent dexterity, yet tough enough to resist punctures and tears. Gloves are 100% nitrile and contain no fillers, silicones, or plasticizers. Ambidextrous.

Approval: CE

Compliance: FDA regulations for food contact



81602-55



### Specifications & Ordering Information


| Catalog number              | Size    | Length | Thickness | Qty/bx | Price/box |
|-----------------------------|---------|--------|-----------|--------|-----------|
| <b>Powdered gloves</b>      |         |        |           |        |           |
| <a href="#">KH-81602-40</a> | Small   | 9½"    | 4 mil     | 100    |           |
| <a href="#">KH-81602-50</a> | Medium  |        |           |        |           |
| <a href="#">KH-81602-60</a> | Large   |        |           |        |           |
| <a href="#">KH-81602-70</a> | X-large |        |           |        |           |
| <b>Powder-free gloves</b>   |         |        |           |        |           |
| <a href="#">KH-81602-45</a> | Small   | 9½"    | 4 mil     | 100    |           |
| <a href="#">KH-81602-55</a> | Medium  |        |           |        |           |
| <a href="#">KH-81602-65</a> | Large   |        |           |        |           |
| <a href="#">KH-81602-75</a> | X-large |        |           |        |           |

## SAFESKIN® Nitrile Exam Gloves

### Purple color easily distinguishes these gloves from latex

These 100% nitrile gloves are soft, latex-free, powder-free, and nonsticky. Low chemical content reduces risk of allergic reactions. Textured fingertips enhance wet or dry grip. Beaded cuff. Ambidextrous.

Compliance: ASTM D 6319, ISO 9002.



86313-21

### Specifications & Ordering Information

| Catalog number              | Size    | Length | Thickness | Qty/bx | Price/box |
|-----------------------------|---------|--------|-----------|--------|-----------|
| <a href="#">KH-86313-20</a> | Small   | 9½"    | 5 mil     | 100    |           |
| <a href="#">KH-86313-21</a> | Medium  |        |           | 100    |           |
| <a href="#">KH-86313-22</a> | Large   |        |           | 100    |           |
| <a href="#">KH-86313-23</a> | X-large |        |           | 90     |           |
| <a href="#">KH-86313-30</a> | Small   | 12"    | 5 mil     | 50     |           |
| <a href="#">KH-86313-31</a> | Medium  |        |           |        |           |
| <a href="#">KH-86313-32</a> | Large   |        |           |        |           |
| <a href="#">KH-86313-33</a> | X-large |        |           |        |           |