

**buy more,
save more**

6-29 packs **save 10%** 30-119 packs **save 15%** 120+ packs **save 20%**

Mix and match packs of any items from the "blue" tables to save!

Bags

BA

WHIRL-PAK® Sterile Sampling Bags, LDPE

Use these sterile LDPE bags to safely hold your solids and liquids. To use, simply flip bag around wire/tape and bend wire to seal. Puncture-proof tabs ensure wire will not poke through even after extensive use.

Sampling Bags are available with or without a white labeling area.

"Speci-Sponge" Bags are designed to sample surfaces and equipment for *Listeria*, *Salmonella*, and other food-borne pathogens; contain a 1½" x 3" sponge that is free of bactericides and noninhibitory; with white labeling area.

Sodium Thiosulfate Bags are EPA approved for potable water collection; contain a nontoxic, nonnutritive 10-mg active ingredient tablet of sodium thiosulfate to neutralize the chlorine at the point of collection; with white labeling area.

Retain Bags are made from a polyethylene/polyester barrier film which has a low permeability for oxygen and moisture; with white labeling area. Designed to safely store your sample for up to two years.

Pocket Bag features an outside pocket (5½" x 5½") to keep your important documentation dry, clean, and separate from the sample.

Black Photo-Sensitive Bag is engineered to keep your sample secure and keep the light out—ideal for pharmaceuticals, food supplements, and other liquids or solids that are sensitive to light. With white labeling strip.

Biohazard Bag features biohazard imprint required by OSHA when storing or transporting hazardous specimens; has unique outside pouch to separate important documents.



Mailing tube
06498-14



Catalog number	Capacity oz (mL)	Size (W x L)	Thickness (mil)	Qty/ pk	Price/ 1 pk	Save 10% 6-29 pks	Save 15% 30-119 pks	Save 20% 120+ pks
Sampling bags without white labeling area								
R-06499-10	4 (118)	3" x 7¼"	2.25	500				
R-06499-20	18 (532)	4½" x 9"	2.5	500				
R-06499-30	24 (710)	6" x 9"	3.0	500				
R-06499-40	36 (1065)	5" x 15"	3.0	500				
R-06499-50	42 (1242)	6" x 15"	3.0	500				
Standup sampling bags with white labeling area								
R-06499-16	4 (118)	3" x 7¼"	2.5	500				
R-06499-21	18 (532)	4½" x 9"	3.0	500				
R-06499-31	24 (710)	6" x 9"	3.0	500				
Sampling bags with white labeling area								
R-06499-60	2 (58)	3" x 5"	2.5	500				
R-06499-70	4 (118)	3" x 7¼"	2.25	500				
R-06499-80	18 (532)	4½" x 9"	2.5	500				
R-06499-85	24 (710)	6" x 9"	3.0	500				
"Speci-Sponge" bags for surface sampling with sponge								
R-06499-13	18 (532)	4½" x 9"	2.5	100				
R-06499-14	55 (1627)	7½" x 11¾"	4.0	100				
R-06499-52†	18 (532)	4" x 9"	2.5	100		—	—	—
Sodium thiosulfate bags for water sampling with tablet								
R-06499-90	3¾ (100)	3" x 7¼"	2.25	100				
R-06499-95	10½ (300)	4½" x 9"	2.5	100				
Retain bags								
R-06498-00	4 (118)	3" x 7¼"	2.5	500				
R-06498-02	18 (532)	4½" x 9"	2.5	500				
Pocket bag								
R-06498-04	24 (710)	6" x 9"	2.5	500				
Black photo-sensitive bag								
R-06498-06	18 (532)	4½" x 9"	2.5	500				
Biohazard bag								
R-06498-08	24 (710)	6" x 9"	2.5	500				

† Model 06499-52 has a sterile glove.

Accessories

R-66203-10 Carrying rack for 2- and 4-oz bags. Double grid design has a second grid approximately 1" below the top to keep bags upright. Rack measures 9¼"L x 6¾"W x 3"H and has 15 compartments that measure 1½"L x 2"W. Vinyl coated, autoclavable.....

R-66203-00 Carrying rack for 18-oz bags. Double grid rack has a second grid approximately 1" below the top to keep bags upright. Rack measures 13¼"L x 10¾"W x 4"H and has 20 compartments that measure 2½"L x 2½"W. Vinyl coated, autoclavable.....



06499-16 06499-21



Water scoop
06498-12



Rack 66203-00

R-06498-12 Water scoop, 18"L coated handle. Two alligator clips hold a variety of bag sizes.....
R-06498-14 Mailing tube for 3¾-oz (100-mL) sodium thiosulfate bag. Polyethylene, 5"L x 2" dia.....