

buy more,
save more

6-29 packs save 10% 30-119 packs save 15% 120+ packs save 20%

Mix and match packs of any items from the "blue" tables to save!



06007-00

Graduated Beaker, Paper

Unwaxed, plain-finish paper beakers have both metric and English graduations (cc and oz). Order dispenser to hold paper beakers for quick and easy access.

Catalog number	Capacity cc (oz)	Qty/ pk	Price/ 1 pk	Save 10% 6-29 pks	Save 15% 30-119 pks	Save 20% 120+ pks
R-06007-00	120 (4)	100				

R-06007-07 Wall-mount dispenser for beaker 06007-00.....



06006-98 06006-92

Stackable Beakers, PP

Polypropylene beakers are graduated in both metric and English (cc and oz) units. Beakers are flexible—a gentle squeeze creates a spout anywhere on the rim. Stacking ribs prevent beakers from sticking together. Withstand temperatures from 0 to 121°C (32 to 250°F).

Autoclavable.

Catalog number	Capacity cc (oz)	Qty/ pk	Price/ 1 pk	Save 10% 6-29 pks	Save 15% 30-119 pks	Save 20% 120+ pks
R-06006-92	120 (4)	1000				
R-06006-98	500 (16)	500				

New



Tri-corner beakers

Cole-Parmer Tri-Corner Beakers, PP

Polypropylene beakers feature a tri-corner design that provides three dripless pouring spouts. Beakers are perfect for use with stirrers—flat bottom keeps stir bars centered. Graduated with an accuracy of ±5%. Withstand temperatures to 121°C (250°F).

Autoclavable.

Catalog number	Capacity mL (oz)	Qty/ pk	Price/ 1 pk	Save 10% 6-29 pks	Save 15% 30-119 pks	Save 20% 120+ pks
R-06010-37	50 (1.5)	100				
R-06010-47	100 (3.5)					
R-06010-57	250 (8)					
R-06010-67	400 (13)	100				
R-06010-77	800 (26)	100				
R-06010-87	1000 (32)	100				
R-06010-05	Five of each above	30				



61003-20 61003-80

Cole-Parmer Economy Beakers, PP

Molded-in graduations are easy-to-read and won't wear off. Polypropylene beakers feature a stable footed base and a no-drip spout. Autoclavable.

Catalog number	Capacity mL (oz)	Qty/ pk	Price/ 1 pk	Save 10% 6-29 pks	Save 15% 30-119 pks	Save 20% 120+ pks
R-61003-12	50 (1.5)	12				
R-61003-22	100 (3.5)	12				
R-61003-42	250 (8)	6				
R-61003-52	400 (13)	6				
R-61003-62	600 (20)	4				
R-61003-72	1000 (32)	3				
R-61003-82	2000 (64)	1				
R-61003-92	4000 (128)	1				



06806-00

06805-00

BenchGuard™ Absorbent Bench Protectors

A highly absorbent upper surface and a tough polycoated base contain spills and protect work surfaces. Choose BenchGuard Extra for twice the absorbency in your most demanding applications. Available in a convenient tear-off roll dispenser or in a box of 50 precut sheets.

Catalog number	Description	Absorbency	Package	Size	Price
ED-06805-00	BenchGuard	400 mL/m ²	Roll	19.3"W x 164 ft L	/rl
ED-06806-00			50 sheets	19.3"W x 2 ft L	/bx
ED-06805-10	BenchGuard Extra	800 mL/m ²	Roll	19.3"W x 164 ft L	/rl
ED-06806-10			50 sheets	19.3"W x 2 ft L	/bx