



Test Tube Brushes

These brushes are specifically designed for cleaning test tubes. Radial-tip brushes have formed brush

sprays at the tip allowing thorough cleaning of all sides. Stainless steel (SS) wire stem radial-tip brushes are extremely durable and provide superior performance. Sponge-tip brushes are made to expand to give a more thorough cleaning while also protecting from scratches. Tied-tip brushes are designed to enter small openings of bottles and tubes, then expand to clean the sides and bottom without scratching. Adjustable brushes spread apart to fit all sizes of test tubes.



Catalog number	Description	Brush bristles			Handle length	Qty/ pk	Price/ pk
		Size	Material	Color			
R-84550-52	Radial tip	4" L x 3/4" W	Nylon	White	6"	12	
R-84551-50		3" L x 1/2" W			5"	12	
R-84551-51		3" L x 1" W			7"	12	
R-84551-52		4" L x 1 3/8" W	Natural	White	9"	12	
R-84551-53		3" L x 5/8" W			9"	12	
R-84551-54		3" L x 7/8" W			9"	12	
R-84551-55	Radial tip/ SS wire stem	4" L x 1 3/8" W	Nylon	White	7"	12	
R-84551-56	3" L x 1/2" W	5"			12		
R-84551-57	3 1/2" L x 3/4" W	5 1/2"			12		
R-84551-58	Sponge tip	4" L x 1 3/8" W	Natural	White	7"	12	
R-84551-59		3" L x 1 3/8" W			6"	12	
R-84551-60		3" L x 1 3/8" W			10"	12	
R-84551-61	Tied tip	3" L x 1 1/2" W	Natural	White	7"	12	
R-84551-62		3 1/2" L x 1/2" W			6"	12	
R-84551-63		6" L x 3/4" W			6"	12	
R-84551-64	Adjustable	3" L x 1 1/4" W	Nylon	White	9"	12	
R-84551-65		6 3/4" L x 1/2" W			2 1/4"	1	

Cylinder/Funnel Brushes

These brushes feature soft wirebound bristles, a tufted end, and smaller diameter plastic handle—ideal for cleaning glass and plastic cylinders and funnels. Handle has hanging hole.



Catalog number	Description	Brush bristles			Handle length	Qty/ pk	Price/ pk
		Size	Material	Color			
R-84550-71	Cylinder/ funnel brush	16" L x 3/8" W	Polyester	White	10 1/2"	12	
R-84550-73		24" L x 3/8" W			17 1/2"	12	

Fine Wire Cleaning Kit

Use this kit to clean syringes, needles, and more. Includes ten 0.15-mm SS wires, ten 0.25-mm SS wires, ten 0.50-mm SS wires, and a pin vise to hold wire during use.



Catalog number	Description	Price
ED-34650-00	Fine wire cleaning kit	

Special-Purpose Brushes


Use these brushes to clean hard-to-reach areas on taps, spouts, plungers, tubes, blades, and equipment fittings. Special design reduces cross-contamination and decreases growth of bacteria and molds by reaching areas prone to growth and trapping particulates.



Catalog number	Description	Brush bristles			Handle length	Qty/ pk	Price/ pk
		Size	Material	Color			
R-84550-57	Tap head	2 1/2" L x 1/4" W	Polyester	White	5 1/2"	12	
R-84550-58	Tap head and tube brush	1/8" tap ; 1/4" tube	Nylon	White	5 1/3"	12	
R-84550-59		1/4" tap ; 1/2" tube			4 1/4"	12	

Burette/Pipette Brushes

Long, flexible handle allows you to easily clean burettes and pipettes.



Catalog number	Description	Brush bristles			Handle length	Qty/ pk	Price/ pk
		Size	Material	Color			
R-84550-56	Burette brush	5" L x 5/8 to 7/8" W	Natural	Black	31"	12	
R-84550-60	Pipette brush	6" L x 3/16" W	Nylon	White	9"	12	
R-84550-61		7 1/2" L x 1/8 to 3/8" W			Nylon	White	9"
R-84550-62		7 1/2" L x 3/16 to 3/4" W	9 1/2"	12			
R-84550-63		7 1/2" L x 3/16 to 5/16" W	9 1/2"	12			
R-84550-64		1 1/2" L x 1/8" W tip ; 4" L x 3/4" W body			12"	12	

Centrifuge and Micro Tube Brushes

Brush design is ideal for cleaning 0.2- to 100-mL centrifuge tubes, vials, bottles, tubes, and a variety of other containers.



Catalog number	Description	Brush bristles			Handle length	Qty/ pk	Price/ pk
		Size	Material	Color			
R-84551-90	Centrifuge tube brush, tied tip	3" L x 3/16-1/2" W	Natural/ nylon	White	6"	12	
R-84551-91	Centrifuge tube brush, radial tip	3 1/2" L x 3/8-1 1/2" W	Natural	White	5 1/2"	12	
R-84551-92	Centrifuge tube brush, radial tip	3" L x 1/2-1 3/16" W			6"	12	
R-84551-93	Centrifuge tube brush, tied tip	4 1/2" L x 3/4-1 1/2" W			7 1/2"	12	
R-84551-94	Micro tube brush	3" L x 1/4" W	Nylon	White	3"	12	
R-84551-95	Micro tube brush, fan tip	3" L x 3/16-1/2" W			4"	1	