


**A**  
00946-03

**B**  
00946-07

**C**  
00946-14

**D**  
00946-18

## Stainless Steel (SS) Laboratory Carts

Stainless steel carts withstand chemical spills and can be steam cleaned

### A Equipment Carts

Models 00946-01 and -03 have a 22-gauge steel top and 20-gauge steel legs; models 00946-04 and -06 have 20-gauge steel tops and 18-gauge steel legs; models 00946-02 and 00946-05 have an 18-gauge steel top and 16-gauge steel legs. Available with or without an equipment drawer.

Catalog number	Shelf area (L x W)	No. of shelves	Shelf clearance	Capacity lb (kg)	Caster dia	Dimensions (L x W x H)	Price
<b>Equipment carts</b>							
R-00946-01				200 (91)	3½"	18¾" x 16¾" x 30"	
R-00946-02	18¾" x 16¾"	3	10¾"	400 (181)	4"	18¾" x 16¾" x 30¾"	
<b>Equipment carts with equipment drawers</b>							
R-00946-03	18¾" x 16¾"	3	8½"	200 (91)	3½"	18¾" x 16¾" x 30"	
R-00946-04		2	17½"			18¾" x 16¾" x 35"	
R-00946-05	18¾" x 16¾"	3	8½"	400 (181)	4"	18¾" x 16¾" x 30¾"	
R-00946-06	36 x 20		10½"			36¾" x 20½" x 35"	

### B General-Purpose Carts

Move laboratory equipment safely and easily. Carts 00946-07 and -09 have 22-gauge steel shelves and 20-gauge steel legs; all other carts have 20-gauge steel shelves and 16-gauge steel legs. Model 00946-12 has two fixed and two swivel casters; all other models feature swivel casters. Available with or without guard rails.

Catalog number	Shelf area (L x W)	No. of shelves	Shelf clearance	Capacity lb (kg)	Caster dia	Dimensions (L x W x H)	Price
<b>General-purpose carts with guard rails</b>							
R-00946-07	24" x 15½"		10"	200 (91)	3½"	27¾" x 16¾" x 33¼"	
R-00946-08	27" x 18"	3	9¾"	400 (181)	4"	31¼" x 19¼" x 33¾"	
<b>General-purpose carts without guard rails</b>							
R-00946-09	24" x 15½"		11¾"	200 (91)	3½"	27¾" x 16¾" x 31½"	
R-00946-10	27" x 18"		11½"	400 (181)	4"	31¼" x 19" x 33¼"	
R-00946-11	35" x 21"	3				39" x 22½" x 37½"	
R-00946-12	49" x 21"		13"	400 (181)	5"	54" x 22½" x 37½"	

### C Heavy-Duty Carts with Guard Rails

Carts feature guard rails to keep loads on shelves and bumpers on legs and handles to protect the surroundings. Shelves are made of 18-gauge steel; 1" square tube legs are made of ½"-thick steel.

Catalog number	Shelf area (L x W)	No. of shelves	Shelf clearance	Capacity lb (kg)	Caster dia	Dimensions (L x W x H)	Price
R-00946-13	27" x 18"	2	17½"			33¾" x 19½" x 34½"	
R-00946-14	21" x 33"	2	18¾"	650 (295)	5"		
R-00946-15		3	9¾"			38¾" x 22¾" x 37¼"	

### D Ergonomic Carts

Designed to reduce body strain. Vertical push handles keep hands in a comfortable position, reducing stress on arms and wrists. Padding on the handles (except models 00946-16 and -17) absorbs vibrations and ensures a good grip. Bottom shelf is placed 20" from the floor to prevent potential strain.

Catalog number	Shelf area (L x W)	No. of shelves	Shelf clearance	Capacity lb (kg)	Caster dia	Dimensions (L x W x H)	Price
R-00946-16	24" x 15½"		12¾"	300 (136)	4"	28" x 16½" x 42¾"	
R-00946-17	27" x 18"					31½" x 19" x 42¾"	
R-00946-18	27" x 18"	2	14"			31½" x 19" x 44½"	
R-00946-19	35" x 21"		16"	500 (227)	5"	39½" x 22" x 44½"	